


SCIENTECH®

ZETA SERIES



Patent
Issued

New Weighing
Technology



New
5-year
Warranty

ISO 9001
REGISTERED

The **MEASURE**
of Confidence

ZETA SERIES

A new, revolutionary weighing technology has recently been developed by Scientech. Since this advanced technology is so exceptional, the US Patent Office has issued Scientech a patent covering this impressive new design. In addition, this innovative technology is now the foundation for a new series of electronic balances which carry Scientech's "Zeta" trademark.



The heart of the new Zeta Series is based on a direct-loading, electromagnetic, force-restoration cell concept, which differs dramatically from normal lever-configuration sensor designs. Thus Zeta models have a more elegant design with fewer moving parts. The result is a more rugged balance with an easily repairable, modular sensor while preserving high performance.

In addition, each unrivaled Zeta model is handcrafted in the USA and retains all the accuracy, durability, and use-friendly standard features you have come to expect from Scientech. The benefits inherent in the Zeta Series also empower Scientech to offer you a new North American 5-year warranty backed by Scientech's ISO9001 certification.




ZETA SERIES STANDARD FEATURES

- ▶ GLP/GMP/ISO compatible
- ▶ Automatic calibration programs
- ▶ RS232 data interface with user selectable baud rates up to 19,200
- ▶ Automatic zero tracking
- ▶ Automatic system checkout sequencing
- ▶ Built-in security bracket
- ▶ Stability indicator
- ▶ Live animal weighing
- ▶ Percent weighing
- ▶ Capacity tracker
- ▶ Microprocessor control
- ▶ Exclusive real time temperature compensation software
- ▶ Parts counting
- ▶ Check weighing
- ▶ Selectable filtering
- ▶ Multiple units of measure
- ▶ Floating Tare™ (dual range models only)
- ▶ Autoranging (dual range models only)


Special FEATURES

Capacity Tracker




This dynamic graphical display automatically depicts the used and unused portions of your weighing range. The tracker is indispensable in taring out container weight and then weighing several components sequentially. You always know how much capacity remains.

Stability Indicator



The balance tells you when the reading is valid, insuring reliable results. Stability sensitivity is user selectable which allows you to customize your balance setup. You may disable the indicator at any time.

Check Weighing




Here is Scientech's answer to a go/no-go weight gauge. The balance tells if your sample weight is too high, too low, or OK. Interpretive error is eliminated. This application is ideal for JIT environments where quality of each sample must be assured.


Floating Tare™ (Dual range models only.)

Dual range balances offer you two balances in one. Used hand-in-hand with autoranging, Floating Tare™ provides low range readability throughout the entire high range capacity of the balance. Access Floating Tare™ with the press of a button and you have ten times more resolution anywhere within the high range domain.


Percent Weighing



By quickly selecting the % mode, you are able to weigh a sample in percent in order to easily determine percent moisture loss, percent solids, or compounding by percent.




Parts Counting



You can even use a Scientech balance for piece counting. Every model empowers you to count very small to large components quickly and accurately. The minimum fixed sample number is ten with a user selectable larger sample size feature at your fingertips.

Vibration Filters



This feature authorizes you to select a filter setting which is best suited to your task. If your surroundings are vibration-free, set the filtering to low which yields results in the fastest possible time. In more adverse conditions, select either the normal or high filter parameter. These selections will insure reliable readings during your different weighing applications.

Units of Measure

These multifunctional balances present a choice of eight different standard units of measure: grams (g), carats (ct), ounces (oz), troy ounces (ozt), pennyweight (dwt), milligrams (mg), kilograms (kg), and pounds (lb). Just push the mode button to select your desired unit of measure. Also available as an option is a custom weight register which allows any unit of measure you wish.

Live Animal Weighing



Weighing using the time averaging mode is ideal for live animal or severe environment weighing. This feature also includes a reset function which allows you to arrive at the final weight reading in seconds rather than minutes with just the push of a button.

SERIES *Specifications*



MODEL	ZSA 80	ZSA 120	ZSA 210	ZSA 210D
WEIGHING MODE	Single	Single	Single	Dual
CAPACITY/RANGE	g 80	120	210	100/200
READABILITY	g 0.0001	0.0001	0.0001	0.0001/0.001
REPEATABILITY (std dev)	g 0.00015	0.00015	0.00015	0.00015/0.0015
LINEARITY	±g 0.0002	0.0002	0.0002	0.0002/0.0005
STABILIZATION TIME	s ~3	~3	~3	~3
PAN SIZE	in Circular 3.5 diameter			
WINDSCREEN	in 7.5W x 8.25D x 9.5H			
BALANCE HOUSING	in 7.5W x 11.25D x 3.25H			



MODEL	ZSP 150	ZSP 250	ZSP 350	ZSP 500	ZSP 404D	ZSP 510D
WEIGHING MODE	Single	Single	Single	Single	Dual	Dual
CAPACITY/RANGE	g 150	250	350	500	40/400	100/500
READABILITY	g 0.001	0.001	0.001	0.001	0.001/0.01	0.001/0.01
REPEATABILITY (std dev)	g 0.0015	0.0015	0.0015	0.0015	0.0015/0.015	0.0015/0.015
LINEARITY	±g 0.0015	0.0015	0.0015	0.0015	0.0015/0.005	0.0015/0.005
STABILIZATION TIME	s ~3	~3	~3	~3	~3	~3
PAN SIZE	in Circular 4.5 diameter					
WINDSCREEN	in Circular 6 diameter x 3.25H					
BALANCE HOUSING	in 7.5W x 11.25D x 3.25H					



Jewelry Scale

MODEL	ZSL 400	ZSL 600	ZSE 250
WEIGHING MODE	Single	Single	Single
CAPACITY/RANGE	g 400	600	250
READABILITY	g 0.01	0.01	0.001
REPEATABILITY (std dev)	g 0.015	0.015	0.0015
LINEARITY	±g 0.015	0.015	0.0015
STABILIZATION TIME	s ~3	~3	~3
PAN SIZE	in Circular 4.5 diameter		
WINDSCREEN	in Optional		7.5W x 8.25D x 9.5H
BALANCE HOUSING	in 7.5W x 11.25D x 3.25H		

Optional ACCESSORIES

GLP/GMP/ISO Compatible

All Scientech balances are GLP/GMP/ISO compatible. Using a computer and one of our small, optional RS232 printers, you may print a hardcopy of either the GLP/GMP Compliance Report and/or a Balance Test for ISO calibration compliance anytime you desire.



Printers



These small thermal printers utilize the standard RS232 data interface which is in each Zeta Series model. Each printer includes interconnecting cable and is ideal for recording your reading for later review.

Custom Weight Register

This factory installed option allows you to weigh in any unit of measure you desire. It is ideal for results in basis weight, any weight per unit area measurement, or even dollars.

Anti-Vibration Table



Ideal for high vibration environments, this small, compact, vibration plate was designed to fit any Scientech balance or scale.

Other accessories include:

Optional remote weighing, NIST traceable calibration weights, spill covers, weighing bowls, rectangular pans, dust covers and software.

Windscreens

Two windscreens are available for Scientech balances:

- ▶ 9.5" high analytical glass wind-screen with 3 sliding glass doors
- ▶ circular glass windscreen with stainless steel cover



Made in the USA

SCIENTECH

1-800-525-0522

5649 ARAPAHOE AVENUE, BOULDER, COLORADO 80303

PHONE: (303) 444-1361 FAX: (303) 444-9229

EMAIL: inst@scientech-inc.com

WEB SITE: <http://www.scientech-inc.com>

An ISO 9001 Registered Company



ISO 9001
REGISTERED

